

# Advanced Spectral Imaging: Homework 2

Yuya Haga

December 23, 2024

## 1 Tasks #1. Colab/Python and spectral files.

### 1.1 Specim IQ

#### 1.1.1 Preview

A spectral image by SpecimIQ was loaded by the Python script shown in Code 4 and converted to RGB preview by Code 5. The gray scale preview and RGB preview are shown in Figure 1.

Gray preview of the image was created by the following code:

```
1 from pathlib import Path
2
3 path = Path('ASI course 2024/group3/Session1/Specim scanner/
4           Color_checker_8_binning/capture/solutions_scan_0110')
5 spectral_image, envi_header = load_spectral_image(path)
6 plt.imshow(spectral_image[:, :, 100], cmap="gray")
```

Code 1: Show gray scale preview

RGB preview of the image was created by the following code:

```
1 rgb_view = reconstruct_rgb_envi(spectral_image, envi_header)
2 plt.imshow(rgb_view)
```

Code 2: Show RGB preview

Functions used in Code 1 and Code 2 are defined in Code 4 and Code 5.





Figure 1: Previews of Specim IQ

### 1.1.2 Plot pair of spectra in one

Figure 2: Plot pair of spectra in one



## 1.2 Spectral Scanner

Spectral image by Spectral Scanner was also processed by the same Python script. The gray scale preview and RGB preview are shown in Figure 3.

The same codes, with different path in Code 1 and Code 2 were used to create the previews.

### 1.2.1 Preview





Figure 3: Previews of Spectral Scanner

### 1.2.2 Plot pair of spectra in one


Figure 4: Plot pair of spectra in one



## 2 Tasks #2. Open ENVI from Japanese spectral camera

A spectral image by Japanese camera was loaded by the Python script shown in Code 6. The gray scale preview and RGB preview are shown in Figure 5.


Figure 5: Gray scale preview of Japanese camera



### 3 Tasks #3. Open ENVI from other byte order

A spectral image of coin was loaded by the Python script shown in Code 7. The gray scale preview and RGB preview are shown in Figure 6. Function `load_envi_header` is defined in Code 4.

Figure 6: Previews of coin



### 4 Tasks #4. Save ENVI spectral image with BIP


Original image was loaded from the ENVI format with the interleave format BIL. The image was saved with the interleave format BIP. The Python script is shown in Code 8.

Function `load_spectral_image` in 8 is defined in Code 4.

The header file of the BIP image is shown in Code 3. There are two changes from original one: description and interleave method.

The preview by FreeLook is shown in Figure 7.

Figure 7: Preview of BIP image by FreeLook



```

1 ENVI
2 description = {Data transformed to BIP}
3 samples = 512
4 lines = 512
5 bands = 204
6 header offset = 0
7 file type = ENVI
8 data type = 12
9 interleave = BIP
10 sensor type = SPECIM IQ
11 byte order = 0
12 default bands = {70,53,19}
13 latitude = 0.00000000
14 longitude = 0.00000000
15 acquisition date = 29-09-2020
16 errors = none
17 binning = {1,1}
18 tint = 121
19 fps = 8.26446
20 wavelength = {
21     397.32,
22     400.20,
23     403.09,
24     405.97,
25     408.85,
26     ...
27 }

```

Code 3: Header file of BIP

## 5 Codes

The following Python scripts were used to complete the tasks. All of those codes are available in the GitHub repository <sup>1</sup>.

```

1 from pathlib import Path
2
3 import numpy as np
4
5
6 def parse_envi_header(lines: list) -> dict[str, str]:
7     """
8     Parses ENVI file content into a structured dictionary

```

<sup>1</sup><https://github.com/8gaU8/ASI-Homeworks>

```

9      This code was written with Github Copilot
10     """
11     envi_data = {}
12     in_block_key = None
13
14     for line_org in lines:
15         line = line_org.strip()
16
17         # Skip empty lines
18         if not line:
19             continue
20
21         # Handle multiline blocks
22         if in_block_key:
23             if line.endswith("}"):
24                 # Closing multiline
25                 envi_data[in_block_key] += line[:-1].strip()
26                 in_block_key = None
27             else:
28                 # Continue multiline
29                 envi_data[in_block_key] += line
30             continue
31
32         # Key-value pair parsing
33         if "=" in line:
34             key, value = map(str.strip, line.split("=", 1))
35             key = key.lower().replace(" ", "_") # Normalize key format
36
37         # Handle block values
38         if value.startswith("{"):
39             # Handle multiline block
40             in_block_key = key
41             # Remove opening '{'
42             value = value[1:].strip()
43             if value.endswith("}"):
44                 # Single-line block
45                 envi_data[key] = value[:-1]
46                 in_block_key = None
47             else:
48                 # Start multiline block
49                 envi_data[key] = value
50         else:
51             # Single-line value
52             envi_data[key] = value
53     return envi_data
54
55
56 def load_envi_header(hdr_file: Path) -> dict[str, str]:
57     """Loads ENVI header file."""
58     with hdr_file.open(encoding="utf-8") as f:
59         header_content = f.readlines()
60         envi_header = parse_envi_header(header_content)
61     return envi_header
62
63
64 def load_spectral_image(file_stem: Path) -> tuple[np.ndarray, dict[str, str]]:
65     """Loads spectral image from ENVI format."""
66     # Load ENVI header

```

```

67     hdr_file = file_stem.with_suffix(".hdr")
68     env_header = load_envi_header(hdr_file)
69
70     # Load parameters
71     interleave = str(env_header["interleave"])
72     lines = int(env_header["lines"])
73     samples = int(env_header["samples"])
74     bands = int(env_header["bands"])
75     data_type = int(env_header.get("data_type", 12))
76
77     # Map ENVI data type to NumPy dtype
78     data_type_map = {
79         1: np.uint8,
80         2: np.int16,
81         3: np.int32,
82         4: np.float32,
83         5: np.float64,
84         6: np.complex64,
85         9: np.complex128,
86         12: np.uint16,
87         13: np.uint32,
88         14: np.int64,
89         15: np.uint64,
90     }
91
92     if data_type not in data_type_map:
93         msg = f"Unsupported data type: {data_type}"
94         raise ValueError(msg)
95
96     dtype = data_type_map[data_type]
97
98     # Load raw data
99     raw_file = file_stem.with_suffix(".raw")
100     with open(raw_file, "rb") as f:
101         raw = np.fromfile(f, dtype=dtype)
102
103     # define shape and transpose order by interleave method
104     if interleave.upper() == "BIL":
105         new_shape = (lines, bands, samples)
106         axis_order = (0, 2, 1)
107     elif interleave.upper() == "BIP":
108         new_shape = (lines, samples, bands)
109         axis_order = (0, 1, 2)
110     elif interleave.upper() == "BSQ":
111         new_shape = (bands, samples, lines)
112         axis_order = (0, 2, 1)
113     else:
114         msg = f"Interleave {interleave} not supported."
115         raise ValueError(msg)
116
117     spectral_image = raw.reshape(new_shape)
118     # change axis order to 'lines, samples, bands'
119     spectral_image = np.transpose(spectral_image, axis_order)
120
121     return spectral_image, env_header

```

Code 4: Load ENVI format images

```

1 import numpy as np
2
3
4 def search_closest_index(wavelengths: list[float], target_wavelength: float) ->
    int:
5     """Finds the index of the closest wavelength to the target wavelength from
        value of envi header."""
6     closest_index = int(np.argmin(np.abs(np.array(wavelengths) -
        target_wavelength)))
7     return closest_index
8
9
10 def reconstruct_rgb_envi(
11     spectral_image: np.ndarray,
12     envi_header: dict[str, str],
13     rgb_wavelengths: tuple[float, float, float] = (632.15, 528.03, 443.56),
14 ):
15     # """Reconstructs RGB image from spectral image with given three wavelengths
        ."""
16     # # Get wavelengths from header
17     wavelengths = [float(data) for data in envi_header["wavelength"].split(",")]
18     rgb_indeces = [search_closest_index(wavelengths, wl) for wl in
        rgb_wavelengths]
19     lines, samples, _bands = spectral_image.shape
20     rgb_view = np.empty((lines, samples, 3))
21     # Reconstruct RGB image
22     for idx, ch in enumerate(rgb_indeces):
23         rgb_view[:, :, idx] = spectral_image[:, :, ch] / np.amax(
24             spectral_image[:, :, ch]
25         )
26     return rgb_view

```

Code 5: Make RGB preview from given specral image

```

1 from pathlib import Path
2
3 import matplotlib.pyplot as plt
4 import numpy as np
5
6 path = Path(
7     "Camera from Japan/colorChecker.nh7"
8 )
9
10 # uint16 float32 #count=spatial_pixels*sample_lines*spectral_bands '>u2' numpy.
    uint16
11 with path.open("rb") as fopen:
12     raw_image = np.fromfile(fopen, dtype=np.uint16)# Calculate number of bands
        from size of image
13 samples = 1280
14 lines = 1024
15
16 size_of_image = raw_image.shape[0]
17 bands = size_of_image // (samples * lines)
18 print(f"Number of bands: {bands}")
19
20 # Reshape the image to 3D array by bands, lines, samples
21 spectral_image = np.reshape(raw_image, (lines, bands, samples))
22 spectral_image = np.transpose(spectral_image, (0, 2, 1))

```



```

23
24 fig, axes = plt.subplots(2, 1, tight_layout=True)
25 axes[0].imshow(spectral_image[:, :, 40], cmap="gray")
26 axes[0].set_title("Preview of spectral image at 40th band")
27
28 axes[1].plot(spectral_image[100, 100, :])
29 axes[1].set_title("Spectra at 100th line and 100th sample")
30 axes[1].set_xlabel("Band")
31 plt.show()

```

Code 6: Load Japanese camera

```

1 from pathlib import Path
2
3 import matplotlib.pyplot as plt
4 import numpy as np
5
6 root = Path("Image of coin (Senop camera)")
7 dat_path = root / "HSI_snapshot__20212101144836.dat"
8 header_path = root / "HSI_snapshot__20212101144836.hdr"
9
10 header = load_envi_header(header_path)
11
12 interleave = header["interleave"]
13 samples = header["samples"]
14 lines = header["lines"]
15 bands = header["bands"]
16 print(interleave)
17
18 with dat_path.open("rb") as fopen:
19     raw_array = np.fromfile(fopen, dtype=">u2")
20
21 # Convert the string values of samples, lines, and bands to integers
22 samples = int(samples)
23 lines = int(lines)
24 bands = int(bands)
25
26 # Reshape the raw_array based on BSQ interleave
27 raw_image = raw_array.reshape((bands, samples, lines))
28 raw_image = np.transpose(raw_image, (1, 2, 0))
29
30
31 fig, axes = plt.subplots(1, 3, tight_layout=True, figsize=(15, 10))
32 axes[0].imshow(raw_image[:, :, 4], cmap="gray")
33 axes[0].set_title("Preview of spectral image at 4th band")
34
35 axes[1].plot(raw_image[100, 100, :])
36 axes[1].set_title("Spectra at 100th line and 100th sample")
37 axes[1].set_xlabel("Band")
38
39 # show RGB preview
40 rgb_preview = raw_image[:, :, [1, 4, 6]].astype(np.float32)
41 rgb_preview /= np.max(rgb_preview, axis=(0, 1))
42 rgb_preview *= 2
43 rgb_preview = np.clip(rgb_preview, 0, 1)
44 axes[2].imshow(rgb_preview)
45 axes[2].set_title("RGB preview")
46

```

```
47 plt.show()
```

Code 7: Load spectral image of coin

```
1 from pathlib import Path
2
3 # load BIL spectral image
4 root = Path("Lectures+Exercises/LectureExercise #5, Freelook and bmp/
   Colorchecker 121 ms (for Freelook demo)/capture/")
5 si_path = root / "Colorchecker"
6 bil_spectral_image, bil_envi_header = load_spectral_image(si_path)
7
8 # reshape BIL to BIP
9 bip_spectral_image = bil_spectral_image.transpose(0, 1, 2)
10 # save BIP spectral image
11 bip_spectral_image.flatten().tofile("colorchecker_bip.raw")
```

Code 8: Load spectral image of coin